

# PADnet™+ Arterial:

the easiest method to test  
for Peripheral Arterial Disease [P.A.D.]

**BioMedix**  
VASCULAR SOLUTIONS

## Office arterial testing

- Training takes less than one day
- Uses automated means to obtain Ankle-Brachial Index/Toe-Brachial Index values and Pulse Volume Recording waveforms

## Addresses a serious health problem

8–12 million Americans have P.A.D. Without proper treatment, 30% are likely to die within five years of P.A.D.-related heart attack or stroke.\*

## For technicians

PADnet+ is extremely easy to use. On-site training takes less than one day, and a nationally recognized certification program is included.

## For primary care doctors

The device sends test results instantly via the web for interpretation by a vascular specialist of your choice.

## For vascular specialists

Segmental pressures (ABI), TBI and PVR tests help identify obstructive disease and determine whether medical or surgical treatment is necessary.

## For patients

This non-invasive test can be completed in 15–20 minutes, during a regular office visit. It saves time and travel to a hospital or specialty lab.

\*Source: P.A.D. Coalition



# What Sets BioMedix PADnet™ + Arterial Testing Apart?

## Meets requirements for claims submission

Segmental pressures (ABI), TBI exams and PVR waveforms meet the requirements for CPT® codes 93922, 93923 and 93924 (optional Post Exercise testing), given approved physical signs and symptoms.\*

## Easy to use interface

The test sequence is pre-set and intuitive, which saves time in testing and study interpretation.

## Rich data for better diagnosis

Patient demographic, history and risk factors are combined with full color reports.

## Study transfer

Two mouse clicks sends the test quickly and securely, ensuring both speed and error-free transmission while maintaining HIPAA compliance.

## Automated technology


Push-button technology eliminates the use of Doppler probes to determine blood pressure and waveforms.

## Testing technology that makes house calls

Because it is internet-ready, specialists can review results at their convenience, from any secure internet connection in the world.

## Coordination improves care

PCPs and DPMs coordinate with vascular specialists to optimize patient care.



**Easy to read reports:** These charts show segmental pressures (ABI) and Pulse Volume Recording waveforms.

## PADnet+ product specifications

### Supplied laptop specifications

- Intel® Core 2 Duo or greater
- Minimum 2.0 GHz CPU speed or greater
- Minimum 3 GB RAM or greater
- DVD ROM drive included
- 56K fax modem
- 10/100 mbps network card
- Internet-ready

### Hardware

- PADnet+ device
- Medical grade cart
- Pressure cuff kit
- Laptop

### User interface

- Microsoft Windows XP, Vista or 7
- PADnet+ test software

\*Current Procedural Terminology (CPT®) is copyright 2007 American Medical Association (AMA). All Rights Reserved. No Fee schedules, basic units, relative values or related listings are included in CPT. The AMA assumes no liability for the data contained herein. CPT is a registered trademark of the American Medical Association. References to CPT codes herein are not intended to convey any endorsement or sponsorship by, or affiliation with, the AMA.



BioMedix Vascular Solutions, Inc. | St. Paul, MN 55110 | Tel (651) 762-4010  
Toll Free (877) 854-0014 | info@BioMedix.com | www.BioMedix.com  
©2010 BioMedix Vascular Solutions, Inc. Saint Paul, MN 55110 | ML-032 Rev 3-10



PADnet is Approved

The American College of  
Foot & Ankle Orthopedics  
& Medicine